

Frequently Asked Questions

- **Can I scan Impressions?** YES as well as a prepped model and multi-die
- **How many microns does it scan to?** 20 Microns
- **What is included?** Scanner, pc, monitor, keyboard, mouse, webcam, and 1 day training
- **Do I have to attend training?** It is highly recommended, but it can be done remotely
- **Is there a buyout at the end of the lease?** Yes \$1
- **What can I scan?** All fixed, orthodontic and any removable that can be hand articulated
- **Is it an open system?** No. This is an affordable and simple scan and send only concept. No design software is necessary, very simple
- **Can I accept files from Dr.'s?** You can accept files from any intra oral scanner and simply send it to us using My Sun Account
- **Is there any design software?** No it is a scan and send only system
- **What is the down payment?** One month down payment of \$199
- **Are there any license fees?** No extra fees at all. No maintenance fee, no digital model fee, free unlimited support
- **What is the turnaround time?** 5 days in lab from when we receive your digital file
- **What is the policy on repairs?** We will send you a loaner while we repair your scanner, just pay shipping costs
- **What is the Warranty?** Warranty covers your first year, for parts and labor.

United States
United Kingdom - Sweden - Netherlands -
Hong Kong - Germany - France



866.561.9777

1800 9th Avenue North
St Petersburg, FL 33713

A professional highly efficient
laboratory is
just a click away !!!



\$199
Per Month


SunScan™

IMPRESSION & MODEL SCANNER

www.sundentallabs.com

FEATURES

SunScanM is an affordable, high quality 3D scanner



BRIDGE WITH CROWN

IMPRESSION

IMPLANT

SunScanM is designed for 3D scanning of an impression and model for manufacturing Crowns, Bridges, Implants, Frames, or articulated Sunflex and Acrylic Partial.



3D White Light Scanner
3D white light scanning technology allows the user

to obtain data that is much clearer than from the conventional laser type.



Twin-camera Technology
SunScanM boasts a function for obtaining a single set of data with two cameras.



Easy Connection
SunScanM has simplified connection between


computer and scanner by accessing them through USB cable only.



Closed System
The basic output format of SunScanM is the STL. SunScanM

solution for labs to scan and send only to Sun Dental manufacturing facility.

Digital Workflow



BENEFITS

- **Fast Turnaround Time**
5 days in lab
- **Eliminate Shipping Charges**
and delayed cases due to broken or lost models
- **Quality Control at the Push of a Button**
Verify the quality of the impression before sending it to the manufacturing facility. Communicate with your clients immediately about defective or poor impression.
- **No second impression**
The model data is stored in the cloud. That means you can reuse the impression file without the need to call your customer for another impression.

All cases are created with an epoxy like digital model that will not abrade or break for a better fit. With the scan only scanner, there is no need to learn the CAD design process.

Grow Your Business by becoming a

FULL SERVICE DIGITAL LAB

<p>SunDenture PREMIUM DENTURES</p> <p>Process \$69 Complete \$119</p>	<p>SunCast PARTIALS</p> <p>Standard Alloy \$69 Vitalium 2000+ \$89</p>	<p>Sunflex PARTIALS</p> <p>Sunflex Process \$85 Sunflex Complete \$119</p>
--	---	---

<p>Suntech PFM CROWNS</p> <p>Non Precious \$45 Semi Precious \$69</p>	<p>Suntech ZIRCONIA</p> <p>Full Zirconia \$55 Layered Zirconia \$69</p>	<p>Suntech COPINGS</p> <p>Zirconia Coping \$29 FCZ Unglazed \$39</p>
--	--	---