



SunCerám[®]
TRANSLUCENT ZIRCONIA



**HIGH Translucent
Full Zirconia**



SunCerám[®]
TRANSLUCENT FULL ZIRCONIA

VS.

Traditional Full Zirconia

1) High Translucency (1.4 times more translucency than traditional zirconia)		1) Not translucent
2) Natural Shade Transition (multilayer)		2) Single shade (monolithic)
3) Anterior and Posterior Crowns and Bridges		3) Only posterior crowns and bridges
4) Flexural Strength (Anterior 615 MPa, Posterior 1200 MPa)		4) Flexural strength 1200 MPa
5) Completely chip resistant		5) Completely chip resistant

Accepting Files from Any Digital Scanner



(866) 561-9777

sundentallabs.com

info@sundentallabs.com